

YHT 2009-04-09

Here is the "regular" protocol we use for differentiating brown preadipocytes

Typical Differentiation Time Course:

Day 1- split into **Differentiation media**

2

3

Day 4 - induce confluent cells (**Induction media**)

5

Day 6- put back into Differentiation media

7

8- aspirate off old media and add fresh Differentiation media

9

Day 10- fully differentiated cells ready to be used for: Oil Red O staining, RNA/DNA extraction, Western analysis etc.

DIFFERENTIATION MEDIA—

DMEM high

10% FBS

Insulin (final concentration = 20 nano M) (stock = 10 mg/ml = 1.7 milli M, add 1.2 micro liter stock in 100 ml medium = 83,333x dilution)

T3 (final concentration = 1 nano M) (stock = 10 micro M, 10,000x)

INDUCTION MEDIA—

DMEM high

10% FBS

Insulin (final concentration = 20 nano M)

T3 (final concentration = 1 nano M)

Indomethacin (final concentration = 0.125 milli M) (stock = 0.125 M, 1000X)

Dexamethasone (final concentration = 5 micro M) (stock = 2 mg/ml = 5 mM, 1000X)

IBMX (final concentration = 0.5 milli M) (stock = 0.25 M, 500X)

[note: in Klein et al, JBC:274, pp34795-34802, 1999, the first paper describing this protocol, the Dex concentration is 0.5 micro M]